



## Opal Mounting Press Recirculating / Chiller Unit



Opal 410 Hot Mounting Press coupled to the dedicated Recirculating / Chiller Unit for the rapid cooling of mounted samples. Ideal for use where no mains water supply is available or a more environmentally friendly cooling water solution is preferred. Product code ATZ 0730007



The self contained chiller system & direct plumbing to the mounting press makes connection simple. The operation of the unit is automatically controlled by the mounting press and cooling times are now reduced to typically a couple of minutes. Capacity of the unit is 6 litres.

MetPrep Technical Bulletin 0049 - 2015. For further information please contact [sales@metprep.co.uk](mailto:sales@metprep.co.uk)

MetPrep Ltd - Curriers Close, Charter Avenue, Coventry, CV4 8AW - 02476 421222