




## Brillant 220 Clamp Vice 40 Options



- Z1304032 Clamp Vice 40 - complete with 45 degree rotating plate  
Z1304042 Small Cross Feed Base - Z axis travel 30 mm - **\*Now discontinued\***

Can be bolted directly to the Brillant 220 T slot bed or used with the Easy Clamp system

- ATZ2236008 Easy Clamp Base  
ATZ220/4 Easy Clamp Adapter Plate for small cross feed vice

The Brillant 220 Clamp Vice 40 complete with the rotatable 45 degree angle plate and the Small Cross Feed Base allows safe and accurate positioning of the sample. Fine control in positioning in the Z axis can be accomplished over a travel distance of 30 mm. The system can be bolted directly the Brillant 220 T slot bed or used in conjunction with the Easy Clamp system for quick changing of accessories.

MetPrep Technical Bulletin 0013 - 2019. For further information please contact [sales@metprep.co.uk](mailto:sales@metprep.co.uk)

MetPrep Ltd - Falklands Close, Charter Avenue, Coventry, CV4 8AU - 02476 421222